




SCPM

Model	HP	KW	DNA/DNM	Suct.Max m	Q.Max L/min	H.Max m	Q. L/min	0	20	40	60	75	80	90	100	105	Weight Kg
SCPM128	0.5	0.37	1" / 1"	9	75	18	H(m)	18	16	12	8	2					9
SCPM146	0.75	0.55	1" / 1"	9	90	25		25	23	19	14	9	7	2			13
SCPM158	1	0.75	1" / 1"	9	100	30		30	28	25	20	14	12	8	2		14
SCPM180	1.5	1.1	1" / 1.25" / 1"	9	105	37		37	34	30	24	18	17	12	7	3	19
SCPM200	2	1.5	1" / 1.25" / 1"	9	105	43		43	40	35	28	21	19	13	7	3	20

Overall and mounting dimension

Model	Dimensions (mm)										Package Dimension		
	DNA	DNM	a	b	c	d	e	f	g	h	A	B	H
SCPM128	1"	1"	263	82	45	175	136	10	88	220	295	195	240
SCPM146	1"	1"	294	82	47	175	142	10	90	224	330	210	265
SCPM158	1"	1"	310	82	43	192	181	10	95	248	330	210	265
SCPM180	1" / 1.25"	1"	345	82	56	215	180	12	112	278	380	255	320
SCPM200	1" / 1.25"	1"	350	82	55	230	200	12	120	300	380	255	320

